

# Intraoperative 3D Detects Malpositioning in **1 OUT OF 5 CASES** Which are Not Visible on 2D Fluoros\*



C-InSight®  
Imagine 3D

C-InSight® is an add-on to your existing C-Arm, converting 2D scans to rotatable 3D images intraoperatively



## Better Clinical outcomes

- Reduces revision cases
- Verification and documentation intraoperatively

## Affordable and Cost Effective

- Increase utilization of all existing C-Arms with minimal investment
- Quickly connects to any C-Arm with a video output
- Reduce the need for post-operative CT scans

## Minimal radiation

- C-Arm radiation exposure is only 5%-10% of the radiation level from a CT scan and dedicated 3D arms.
- Reduces radiation exposure according to FDA guidelines for both patients and OR staff




## Simple and Quick

- High Image Quality
- Simple - Plug and Play
- Quick - 2D to 3D under 2 Minutes
- Intuitive - review in any angle and slice by slice

## Converts 2D to 3D Intraoperatively in 3 easy steps:

1. Place the C-InSight sterile target on the patient at the region of interest
2. Acquire a 15-20 second continuous fluoroscopic scan of the patient's anatomy
3. The C-InSight converts the scans within 2 minutes

# CLINICAL VALUE OF 3D IMAGING



▶ Standard 2D C-Arm with C-InSight 3D  
Slice 83 out of 115 slices



▶ Standard 2D C-Arm  
Lateral image

## Intra-articular fracture of lateral epicondyle of humerus

3D imaging detects malpositioning which cannot be viewed with regular 2D imaging. Reviewing from all angles and slice by slice, the surgeon can verify the positioning and correct it intraoperatively.



[www.mazorst.com](http://www.mazorst.com)

**Company Headquarters**  
**Mazor Surgical Technologies Ltd.**

Tel: + 972-(0)4-627-0171  
Fax: +972-(0)4-637-7234  
PO Box 3104, Southern  
Caesarea Park 38900, Israel  
E-mail: [info@mazorst.com](mailto:info@mazorst.com)

**USA Office**  
**Mazor Surgical Technologies Inc.**

4361 Shackleford Road  
Norcross, GA 30093  
Tel: (770)-564-5790  
Fax: (770)-564-5791  
E-mail: [usa@mazorst.com](mailto:usa@mazorst.com)

**European Office**  
**Mazor Surgical Technologies UG**

Borkstr. 10,  
48163 Muenster, Germany  
Tel: +49 251 32266-0  
Fax: +49 251 32266-52  
E-mail: [germany@mazorst.com](mailto:germany@mazorst.com)